## cytoactiv<sup>®</sup> HPV L1 High Risk Set

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
### Literatur

1. McMurray, HR.; Nguyen, D.; Westbrook, TF.; Mcance, DJ.: Biology of Human Papillomavirus. Int. J. Pathol. 2001, 82: 13-33
2. Prömer et al., 2004: Nachweis produktiver HPV-Infektionen am zytologischen Routine-Abstrichpräparat. J. Fertil. Reprod. 4/2004.
3. Melsheimer et al., 2003: Immunocytochemical detection of HPV High risk type L1 capsid protein in LSIL and HSIL as compared with detection of HPV L1 DNA. Acta cytologica, 47: 124-128.
4. Birner et al., 2001: Signal-amplified colorimetric in situ hybridisation for assessment of human papillomavirus infection in cervical lesions. Mod. Pathol. 14(7):702-9.
5. Griesser et al.; 2004: Correlation of Immunochemical detection of HPV L1 capsid protein in Pap smears with regression of High risk positive mild/moderate dysplasia. AQCH, Vol. 26, No. 5: 241-245
6. Griesser et al.: 2006, Immunochemical detection of HPV L1 capsid protein: a prognostic marker for early squamous lesions of the cervix uteri. Abstract Eurogin2006, Paris
7. Griesser, H.; Sander, H.; Hilfrich, R.: HPV L1 capsid protein detection as prognostic marker for hr HPV associated early dysplastic lesions. Pathologica 2006, 98: 368
8. zur Hausen, H.: Papillomaviruses and cancer: from basic studies to clinical application. Nat. Rev. Cancer 2002, 2: 342-350
9. Munoz, N.; Bosch, FX.; de Sanjose, S. et al.: Epidemiologic classification of human papillomavirus types associated with cervical cancer. Engl. J. Med. 2003, 348: 518-527
10. Hilfrich R, Hariri J: Prognostic relevance of HPV L1 capsid protein detection within mild to moderate dysplastic lesions of the cervix uteri in combination with a second biomarker p16. AQCH, 2008, Vol 30, No 2, 78-82
11. Rauber, D.; Mehlhorn, G.; Fasching, P. A.; Beckmann, M. W.; Ackermann, S.: Prognostic significance of the detection of the human papillomavirus L1 protein in smears of mild to moderate cervical intraepithelial lesions. European Journal of Obstetrics and Gynecology and Reproductive Biology, October 2008, Vol.140, Issue 2, pages 258-262
12. Griesser H., Sander H., Walczak C., Hilfrich R.: Immunochemical detection of HPV L1 capsid protein: prognostic marker for early squamous lesions of the cervix, Acta cytologica 2007, Vol 52, 268

Cytoimmun Diagnostics GmbH  
Delaware Avenue 1-3  
D-66953 Pirmasens

Phone 0049 6331 145 333-0  
Fax 0049 6331 145 333-99  
E-mail info@cytoimmun.de